

# Enthalpy and Entropy of Zinc with Copper Sulfate

## INSTRUCTOR RESOURCES

### The CCLI Initiative

#### Learning Objectives

- determine changes in enthalpy and entropy of the reaction of zinc with copper sulfate using two methods: electrochemistry and calorimetry.
- compare the enthalpy values obtained by the two methods.

#### Procedure Overview

- a simple electrochemical cell  $\text{Cu(s)}/\text{CuSO}_4(\text{aq}) // \text{ZnSO}_4(\text{aq}) / \text{Zn(s)}$  is constructed in a Chem-Carrou-Cell™ and voltages are measured at different temperatures.
- a spreadsheet is used to plot a linear regression graph of voltage *versus* temperature, and the graph is used to calculate the enthalpy and entropy changes for the cell reaction.
- heat of reaction of zinc powder with 0.5 M  $\text{CuSO}_4$  solution is determined calorimetrically.

Name \_\_\_\_\_ Section \_\_\_\_\_ Date \_\_\_\_\_

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Report Sheet

#### Part I: Electrochemical Procedure

1. From your graph of Voltage *versus* T, copy the straight line equation you have obtained.
2. From the slope of your equation, calculate  $\Delta S$  for the cell reaction.
3. From the y-intercept, calculate  $\Delta H$  for the cell reaction.
4. Using the entropy and enthalpy changes determined above, calculate  $\Delta G$  at 25 C.
5. Calculate  $\Delta G$  using Equation 2.
6. From your straight line equation, calculate the voltage,  $\mathcal{E}^\circ$ , at 25 °C.

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Report Sheet (page 2)

#### Part I : Calorimetric Method

weight of beaker + solution	_____
weight of empty beaker	_____
mass of the solution	_____
weight of boat + zinc	_____
weight of empty boat	_____
mass of zinc	_____
initial temperature	_____
final temperature	_____

1. Write the equation for the reaction that took place in the calorimeter.
  
  
  
  
  
  
  
  
  
  
2. Calculate the heat of reaction per mole of zinc.
  
  
  
  
  
  
  
  
  
  
3. Compare the values of  $\Delta H$  obtained by the two methods.

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Questions/Problems

1. Compare the calculated value of  $\mathcal{E}$  at 25 °C to the value obtained from your graph at the same temperature.
  - (a) Is it close to the theoretical value of 1.1 V?
  - (b) What might be the reasons for any difference between the theoretical and experimental values of you obtained?
2. A student performs this experiment incorrectly, leaving the alligator clips hooked to the cell for the entire time. Assume that the concentration of  $\text{Cu}^{2+}$  changes from 0.50 M to 0.40 M and that the  $\text{Zn}^{2+}$  concentration changes from 0.50 M to 0.60 M. How much will these concentration changes affect the measured potential?
3. Why is it necessary in this experiment to assume that  $\Delta H$  and  $\Delta S$  for the reaction are independent of temperature?

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Suggested Answers to Questions/Problems

1. Compare the calculated value of  $\mathcal{E}$  at 25 C to the value obtained from your graph at the same temperature.

*Answers will vary.*

- (a) Is it close to the theoretical value of 1.1 V?

*Answers will vary.*

- (b) What might be the reasons for any difference between the theoretical and experimental values of you obtained?

*The original concentrations of the solutions might be incorrect or the alligator clips might be left in place too long causing significant current flow in the solutions.*

2. A student performs this experiment incorrectly, leaving the alligator clips hooked to the cell for the entire time. Assume that the concentration of  $\text{Cu}^{2+}$  changes from 0.50 M to 0.40 M and that the  $\text{Zn}^{2+}$  concentration changes from 0.50 M to 0.60 M. How much will these concentration changes affect the measured potential?

*From the Nernst equation (at 25 )*

$$= - (0.0592/n)\log Q$$

*the concentration term is*

$$(0.0592/2) \log[\text{Zn}^{2+}] / [\text{Cu}^{2+}]$$

*This term is zero for  $[\text{Zn}^{2+}] = [\text{Cu}^{2+}] = 0.50 \text{ M}$ . If the concentrations of  $\text{Cu}^{2+}$  and  $\text{Zn}^{2+}$  change to 0.40 M and 0.60 M, respectively, the potential changes by*

$$(0.0592/2) \log(0.60/0.50) = 0.005 \text{ V}$$

3. Why is it necessary in this experiment to assume that  $\Delta H$  and  $\Delta S$  for the reaction are independent of temperature?

*In order to use the equation*

$$= \frac{\Delta S (T)}{n \mathcal{F}} - \frac{\Delta H}{n \mathcal{F}}$$

*to obtain values of  $\Delta S$  and  $\Delta H$ , we plot  $\mathcal{E}$  versus  $\mathcal{T}$ . In doing so, we assume that  $T$  is the only variable on the right side of the equation.*

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Sample *MicroLAB* Program

#### Part I: Electrochemistry

Experiment name: *delta G experiment*

Sensors: **Temperature probe**: X axis, Col. A, DD on top, units = °C; **Voltage**: Y1 axis, Col B, DD on bottom, units = volts.

Special program:

Read Sensors

Repeat every 0.5 seconds

    If Delta TempIC > +/- 5.00

*Read Sensors*

*Else*

*End If*

Until **Stop** Button is pressed

#### Part II: Calorimetry

Experiment file name: *heat of solution*

Sensors used: **Timer 1**: X axis, Col. A, DD on top, units = sec; **Temperature probe**: Y axis, Col. B, DD on bottom, units = °C;

Program

*Repeat every 0.5 seconds*

*Read Sensors*

    Until **Stop** button is clicked

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Tips and Traps

1. For the electrochemical section of the experiment, the best results were obtained below 40 °C.
2.  $\Delta S$  values are quite inconsistent. They depend on the slope reading of the voltage vs. a temperature graph. The slope is always a very small number; therefore, any variation in slope from one trial to another would change the entropy value greatly.
3. The experiment can be performed in a Styrofoam cup. The estimated heat capacity is 10 J °C<sup>-1</sup>.
4. Any combination of small containers can be used in place of the Chem Carron-Cell™.
5. Chem Carron-Cell™ is available from Freeman, Cooper & Co., San Francisco, CA 94133.

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Sample Data

#### Part I: Electrochemical Procedure

1. From your graph of Voltage versus T, copy the straight line equation you obtained

$$V = 0.0003648 (T) + 0.9502$$

2. From the slope of your equation, calculate S for the cell reaction.

$$\Delta S = (2 \text{ mol e}) (96,500 \text{ C/mol e}) (0.0003648 \text{ V/K})$$

$$\Delta S = 70 \text{ J/K}$$

3. From the y-intercept, calculate H for the cell reaction.

$$\Delta H = (2 \text{ mol e}) (96,500 \text{ C/mol e}) (0.9502 \text{ V})$$

$$\Delta H = 1.83 \times 10^5 \text{ J or } 183 \text{ kJ}$$

4. Using of entropy and enthalpy changes determined above, calculate G at 25 C.

$$\Delta G = (183 \text{ kJ}) (298 \text{ K}) (0.070 \text{ kJ/K})$$

$$\Delta G = 204 \text{ kJ}$$

5. Calculate  $\Delta G$  using Equation 2.

$$\Delta G = (2 \text{ mol e}) (96,500 \text{ C/mol e}) (1.06 \text{ V})$$

$$\Delta G = 2.05 \times 10^5 \text{ J or } 205 \text{ kJ}$$

6. From your straight line equation, calculate the voltage,  $\mathcal{E}$ , at 25 C.

$$= 0.0003648 (298) + 0.9502$$

$$= 1.06 \text{ V}$$




## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Sample Data (page 2)

#### Part II: Calorimetric Method

weight of beaker + solution	<u>119.8</u> g
weight of empty beaker	<u>65.6</u> g
mass of the solution	<u>54.2</u> g
weight of boat + zinc	<u>0.8975</u> g
weight of empty boat	<u>0.4043</u> g
mass of zinc	<u>0.4932</u> g
initial temperature	<u>24 C</u>
final temperature	<u>30 C</u>

1. Write the equation for the reaction that took place in the calorimeter.



2. Calculate the heat of reaction per mole of zinc.

$$\Delta H = (3.8 \text{ J g}^{-1} \text{ C}^{-1} \times 54.2 \text{ g} + 30 \text{ J C}^{-1}) (24 \text{ C} - 30 \text{ C}) / (0.4932 \text{ g} / 65.38 \text{ g/mol})$$

$$\Delta H = 1.88 \times 10^5 \text{ J/mol or } 188 \text{ kJ}$$

3. Compare the values of  $\Delta H$  obtained by the two methods.

*The electrochemical value of 183 kJ compares very well with the calorimetric value of 188 kJ.*

## ENTHALPY AND ENTROPY OF ZINC WITH COPPER SULFATE

### Laboratory Preparation (per student station)

#### Part I: Electrochemistry

##### Equipment

- 600 ml beaker
- 50 ml beaker
- tweezers
- Chem-Carron-Cell™
- one pair of black and red alligator clip leads

##### Supplies

- filter paper strips (0.5 x 5 cm)
- sandpaper
- ceramic tile

##### Chemicals

Exact quantities needed are listed below. A minimum 50% excess is recommended.

- copper metal strip (0.5 x 5 cm)
- zinc metal strip (0.5 x 5 cm)
- 0.5 M CuSO<sub>4</sub> (5 ml)
- 0.5 M ZnSO<sub>4</sub> (5 ml)
- 0.1 M KNO<sub>3</sub> (5 ml)

#### Part II: Calorimetry

##### Equipment

- thermistor
- 600 ml beaker
- 400 ml beaker
- 150 ml beaker
- 50 ml graduated cylinder
- insulated cover (for 100 ml beaker)

##### Supplies

- towel

##### Chemicals

Exact quantities needed are listed below. A minimum 50% excess is recommended.

- 0.5 M CuSO<sub>4</sub> (50 ml)
- zinc powder (0.5 g)

##### Safety and Disposal

dispose of the wastes into specially marked containers.