

## Hot Packs, Cold Packs, and Heats of Solution

### INSTRUCTOR RESOURCES

#### The CCLI Initiative

#### Learning Objectives

- use the *MicroLAB* interface to measure the heats of solution ( $\Delta H_{\text{soln}}$ , in kJ/mol) of several salts.
- calculate the heats of formation ( $\Delta H_f$ , in kJ/mol) of single aqueous ions.
- predict the heats of solution of additional salts.
- design a hot pack and a cold pack for specified temperature changes.

#### Procedure Overview

- the *MicroLAB* interface is used to measure temperature changes in a Styrofoam cup calorimeter.
- $\Delta T$  values are determined graphically for the dissolution of four salts.
- $\Delta H_{\text{soln}}$  calculated for each salt.
- $\Delta H_f$  calculated for the ions in each of the four salts.
- $\Delta H_f$  values for the ions used to find  $\Delta H_{\text{soln}}$  for two additional salts.
- a cold pack designed using  $\text{NH}_4\text{NO}_3$ , and a hot pack designed using  $\text{MgSO}_4$ .

**HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION****Report Sheet****Measurement of  $\Delta H_{\text{soln}}$  (data and calculations):**

Salt	Mass (g)	$T_i$ (°C)	$T_f$ (°C)	$\Delta T$ (°C)	$\Delta H_{\text{soln}}$ (kJ/mol)
MgSO <sub>4</sub>	_____	_____	_____	_____	_____
K <sub>2</sub> SO <sub>4</sub>	_____	_____	_____	_____	_____
KNO <sub>3</sub>	_____	_____	_____	_____	_____
NH <sub>4</sub> NO <sub>3</sub>	_____	_____	_____	_____	_____

Sample calculation for:

 $\Delta H_{\text{soln}}$  for MgSO<sub>4</sub> (kJ/mol) $\Delta H_{\text{soln}}$  for K<sub>2</sub>SO<sub>4</sub> (kJ/mol) $\Delta H_{\text{soln}}$  for KNO<sub>3</sub> (kJ/mol) $\Delta H_{\text{soln}}$  for NH<sub>4</sub>NO<sub>3</sub> (kJ/mol)

**HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION****Report Sheet (page 2)****Calculated  $\Delta H_f$  of ions and  $\Delta H_{\text{soln}}$  of salts:**

Ions	$\Delta H_f$ (kJ/mol)	Salts	$\Delta H_{\text{soln}}$ (kJ/mol)
$K^+$	- 251.2	$Mg(NO_3)_2$	_____
$Mg^{2+}$	_____	$(NH_4)_2SO_4$	_____
$NH_4^+$	_____		
$SO_4^{2-}$	_____		
$NO_3^-$	_____		

**Calculations for** $\Delta H_f$  for  $Mg^{2+}$ : $\Delta H_f$  for  $NH_4^+$ : $\Delta H_f$  for  $(SO_4^{2-})$  $\Delta H_f$  for  $NO_3^-$  $\Delta H_f$  for  $Mg(NO_3)_2$  $\Delta H_f$  for  $(NH_4)_2SO_4$

**HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION**

**Report Sheet (Page 3)**

**Design of cold and hot packs:**

**Cold pack:**

Calc. grams of  $\text{NH}_4\text{NO}_3$ :

Calc. grams of water:

**Hot Pack:**

Calculate grams of  $\text{MgSO}_4$ :

Calculate grams of water:

<b>Cold Pack</b>	<b>Hot Pack</b>
_____ g $\text{NH}_4\text{NO}_3$	_____ g $\text{MgSO}_4$
_____ g water	_____ g water



## HOT PACKS, COLD PACKS AND HEATS OF SOLUTION

### Tips and Traps

1. Stirring is critical in this lab.
2. Stirring should be vigorous without splashing.
3. Magnetic stirrers may cause interference with the computers. If so, manual stirring will be necessary.
4. Salts must be thoroughly dried before lab, and kept in tight containers to prevent uptake of additional moisture.
5. Make sure that the salts are added to the water in such a manner that they do not stick to the sides of the Styrofoam cup.
6. The thermistor must be *in* the solution without touching the bottom of the container or the magnetic stir bar.
7. An alternative method for supporting the foam cups is a suitable iron ring (3-inch).
8. If temperature drift is observed after salt addition, extrapolation of temperature to the time of mixing will provide a more reliable  $T_f$  value.
9. You need a cover for the foam cups. There are many alternatives; two are suggested below:
  - a. a plastic coffee cup lid.
  - b. a foam square, 3.5" x 3.5" x 0.5" with a hole in the middle for the thermistor. A rubber band or small piece of rubber hose around the thermistor will prevent it from slipping through the foam block or the coffee cup cover.
10. The heat capacity for the calorimeter can be experimentally determined if so desired. A common method is based on the heat produced on mixing known amounts of HCl and NaOH.
11. In designing the hot pack and cold pack, the assumption of  $4.18 \text{ J/g}^\circ\text{C}$  for the *solution* introduces some error because the solution has a different specific heat than pure water.

## HOT PACKS, COLD PACKS AND HEATS OF SOLUTION

### Sample Data (Page 1)

Measurement of  $\Delta H_{\text{soln}}$  (data and calculations):

Salt	Mass (g)	$T_i$ (°C)	$T_f$ (°C)	$\Delta T$ (°C)	$\Delta H_{\text{soln}}$ (kJ/mol)
MgSO <sub>4</sub>	5.00	20.1	27.5	7.4	-82
K <sub>2</sub> SO <sub>4</sub>	5.00	19.8	18.1	-1.7	27
KNO <sub>3</sub>	5.00	20.0	16.2	-3.8	36
NH <sub>4</sub> NO <sub>3</sub>	5.00	20.0	16.4	-3.6	27

Calculations for each of the four salts:

$$\Delta H_{\text{soln}} (\text{MgSO}_4) = - \frac{(7.4^\circ\text{C}) (120 \text{ g/mol}) (463 \text{ J/}^\circ\text{C})}{(5.00 \text{ g}) (1000 \text{ J/kJ})}$$

$$= -82 \text{ kJ/mol}$$

$$\Delta H_{\text{soln}} (\text{K}_2\text{SO}_4) = - \frac{(-1.7^\circ\text{C}) (174 \text{ g/mol}) (463 \text{ J/}^\circ\text{C})}{(5.00 \text{ g}) (1000 \text{ J/kJ})}$$

$$= 27.4 \text{ kJ/mol}$$

$$\Delta H_{\text{soln}} (\text{KNO}_3) = - \frac{(-3.8^\circ\text{C}) (101 \text{ g/mol}) (463 \text{ J/}^\circ\text{C})}{(5.00 \text{ g}) (1000 \text{ J/kJ})}$$

$$= 36 \text{ kJ/mol}$$

$$\Delta H_{\text{soln}} (\text{NH}_4\text{NO}_3) = - \frac{(-3.6^\circ\text{C}) (80 \text{ g/mol}) (463 \text{ J/}^\circ\text{C})}{(5.00 \text{ g}) (1000 \text{ J/kJ})}$$

$$= 27 \text{ kJ/mol}$$

## HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION

### Sample Data (page 2)

Calculated  $\Delta H_f$  of ions and  $\Delta H_{\text{soln}}$  of salts:

Ions	$\Delta H_f$ (kJ/mol)	Salts	$\Delta H_{\text{soln}}$ (kJ/mol)
$K^+$	-251.2	$Mg(NO_3)_2$	<u>-86</u>
$Mg^{2+}$	<u>-464</u>	$(NH_4)_2SO_4$	<u>10</u>
$NH_4^+$	<u>-133</u>		
$SO_4^{2-}$	<u>-904</u>		
$NO_3^-$	<u>-206</u>		

Sample calculations: (a)  $\Delta H_f$  for  $Mg^{2+}$ , (b)  $\Delta H_{\text{soln}}$  for  $Mg(NO_3)_2$

(a) First calculate  $\Delta H_f$  ( $SO_4^{2-}$ )

$$\begin{aligned} \Delta H_f (SO_4^{2-}) &= \Delta H_{\text{soln}} (K_2SO_4) + \Delta H_f^\circ (K_2SO_4) - 2 \Delta H_f (K^+) \\ &= [(27) + (-1433) - 2(-251)] \text{ kJ/mol} \\ &= -904 \text{ kJ/mol} \end{aligned}$$

Next calculate  $\Delta H_f$  ( $Mg^{2+}$ )

$$\begin{aligned} \Delta H_f (Mg^{2+}) &= \Delta H_{\text{soln}} (MgSO_4) + \Delta H_f^\circ (MgSO_4) - \Delta H_f (SO_4^{2-}) \\ &= [(-82) + (-1285) - (-904)] \text{ kJ/mol} \\ &= -464 \text{ kJ/mol} \end{aligned}$$

(b)  $\Delta H_{\text{soln}} (Mg(NO_3)_2) = \Delta H_f (Mg^{2+}) + 2 \Delta H_f (NO_3^-) - \Delta H_f^\circ (Mg(NO_3)_2)$

$$\begin{aligned} &= [(-464) + 2(-206) - (-790)] \text{ kJ/mol} \\ &= -86 \text{ kJ/mol} \end{aligned}$$



## HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION

### Sample Data (page 2 continued)

Calculate  $\Delta H_f(\text{NO}_3^-)$

$$\begin{aligned} &= \Delta H_{\text{soln}}(\text{KNO}_3) + \Delta H_f^\circ(\text{KNO}_3) - \Delta H_f(\text{K}^+) \\ &= [(36) + (-492.7) - (-251) \text{ kJ/mol}] \\ &= -206 \text{ kJ/mol} \end{aligned}$$



Calculate  $\Delta H_f(\text{NH}_4^+)$

$$\begin{aligned} &= \Delta H_{\text{soln}}(\text{NH}_4\text{NO}_3) + \Delta H_f^\circ(\text{NH}_4\text{NO}_3) - \Delta H_f(\text{NO}_3^-) \\ &= [(27) + (-365.6) - (-206) \text{ kJ/mol}] \\ &= -133 \text{ kJ/mol} \end{aligned}$$

Calculate  $\Delta H_{\text{soln}}((\text{NH}_4)_2\text{SO}_4)$

$$\begin{aligned} &= 2 \times \Delta H_{\text{soln}}(\text{NH}_4^+) + \Delta H_f^\circ(\text{SO}_4^{2-}) - \Delta H_f((\text{NH}_4)_2\text{SO}_4) \\ &= [2 \times (-133) + (-904) - (-1179.5) \text{ kJ/mol}] \\ &= 10 \text{ kJ/mol} \end{aligned}$$

Sample *MicroLAB* heating and cooling screens



## HOT PACKS, COLD PACKS, AND HEATS OF SOLUTION

### Sample Data (page 3)

Design of cold and hot packs:

Cold Pack	Hot Pack
$\frac{65}{335}$ g NH <sub>4</sub> NO <sub>3</sub> g water	$\frac{90}{310}$ g MgSO <sub>4</sub> g water

Cold pack:

$$\text{(NH}_4\text{NO}_3) \quad \frac{(4.18 \text{ J/g}^\circ\text{C}) (400 \text{ g}) (13^\circ\text{C}) (80 \text{ g/mol})}{(27 \text{ kJ/mol}) (1000 \text{ J/kJ})}$$

$$= 65 \text{ g NH}_4\text{NO}_3$$

$$\text{(water)} \quad 400 - 65 = 335 \text{ g H}_2\text{O}$$

Hot pack:

$$\text{(MgSO}_4) \quad \frac{(4.18 \text{ J/g}^\circ\text{C}) (400 \text{ g}) (37^\circ\text{C}) (120 \text{ g/mol})}{(82 \text{ kJ/mol}) (1000 \text{ J/kJ})}$$

$$= 90 \text{ g MgSO}_4$$

$$\text{(water)} \quad 400 - 90 = 310 \text{ g H}_2\text{O}$$

## COLD PACKS, HOT PACKS, AND HEAT OF SOLUTION

### Suggested Answers to Questions/Problems

1. Look up the literature value for the heat of solution of each of the four salts that you studied. Determine the percent error between your value and the literature value. Provide some possible reasons for any differences that are observed.

#### *Heats of Solution (from Lange - Handbook of Chemistry)*

<i>Compound</i>	$\Delta H_{soln}$	<i>Your Value</i>	<i>Difference</i>	<i>% error</i>
MgSO <sub>4</sub>	-84.9 kJ/mol			
K <sub>2</sub> O <sub>4</sub>	+26.7 kJ/mol			
KNO <sub>3</sub>	+35.6 kJ/mol			
NH <sub>4</sub> NO <sub>3</sub>	+26.4 kJ/mol			

1. *Heat of solution varies with extent of dilution.*
  2. *The number of significant figures on measured temperature changes as a function of the care in probe calibration.*
  3. *Calorimeter heat capacity used is an average value, which will vary somewhat with each salt.*
2. Obtain information on costs of the chemicals used in the hot pack and cold pack that you designed and determine the total materials cost to produce these items. Compare these costs to the actual costs of hot packs and cold packs available in local stores.

#### *Approximate material costs:*

*Cold Pack - \$0.30 (based on ammonium nitrate at \$8/lb)*

*Hot Pack - \$0.38 (based on magnesium sulfate at \$10/lb)*

#### *Actual costs:*

*Check local sporting goods stores.*

## HOT PACKS, COLD PACKS AND HEATS OF SOLUTION

### Laboratory Preparation (per student station)

#### Equipment

- temperature probe
- ring stand
- utility clamp
- two foam cups (nested into a 400 ml beaker to form a calorimeter)
- cover for foam cup
- magnetic stirrer and stir bar
- 100 ml graduated cylinder
- MicroLab interface
- ice (for calibration)
- hot water (for calibration)

#### Supplies

- weighing paper or weighing boat

#### Chemicals

Each of these salts should be dried. On the morning of the lab, heat them in the oven in a beaker covered with a watch glass for 2 hours at the temperature indicated. Stir frequently to break up the material and obtain complete removal of unwanted water. Keep in tightly capped containers thereafter. Make a set of separate containers for each lab section.

Per pair of students

- 12 g  $\text{MgSO}_{4(s)}$  (-6  $\text{H}_2\text{O}$  @ 150 °C, mp 1124 °C, heat at 155 °C)
- 12 g  $\text{K}_2\text{SO}_{4(s)}$  (mp 1069 °C, heat at 110 °C)
- 12 g  $\text{KNO}_{3(s)}$  (mp 334 °C, heat at 110 °C)
- 12 g  $\text{NH}_4\text{NO}_{3(s)}$  (mp 169.6 °C, heat at 110 °C)

#### Safety and Disposal

- make sure students wear eye protection at all times.
- all resulting solutions may be flushed down the drain with plenty of water.