

Crystal Violet Colorimetry:  
a Beer's Law Investigation

## INSTRUCTOR RESOURCES

The CCLI Initiative

Computers in Chemistry Laboratory Instruction

### Learning Objectives

- illustrate the basic principles of colorimetry.
- demonstrate the components of a colorimeter and how a colorimeter is interfaced to a computer.
- discover the Beer's Law relationship and apply it to the analysis of an unknown solution.

### Procedure Overview

- "ideal" simulated data is manually entered into the *MicroLAB* spreadsheet and a series of plots using different functions operating on the current readings are constructed in order to discover a linear relationship between the current and concentration variables.
- a series of standard crystal violet solutions are prepared and analyzed using the *MicroLAB* colorimeter. A Beer's Law plot of these concentrations is constructed using the relationship discovered with the simulated data.
- an unknown solution of crystal violet is analyzed in the colorimeter and the Beer's Law plot is used to determine its concentration.

**COLORIMETRY: A BEER'S LAW INVESTIGATION****Report Sheet**

Dilution Table for Standards

	#1	#2	#3	#4	#5
<b>Desired CV Molarity</b>	$2.00 \times 10^{-6}$	$4.00 \times 10^{-6}$	$6.00 \times 10^{-6}$	$8.00 \times 10^{-6}$	$10.0 \times 10^{-6}$
<b>mL CV</b>					
<b>mL H<sub>2</sub>O</b>					
<b>Total mL</b>	15.00	15.00	15.00	15.00	15.00

1. Record the number of your unknown \_\_\_\_\_

\_\_\_\_\_ Record the absorbance of your unknown \_\_\_\_\_

\_\_\_\_\_ Record the concentration of your unknown \_\_\_\_\_

2. Show your calculations for the concentration of your unknown.



**COLORIMETRY: A BEER'S LAW INVESTIGATION****Questions/Problems**

1. When iron (III) chloride is added to a solution of salicylic acid a purple solution forms. The photocell current read with the colorimeter for a series of standard solutions is given below. Convert these values to absorbencies. (The blank current is 100.0.)

FeCl <sub>3</sub>	Transmittance	Absorbance
$40.0 \times 10^{-5}$ M	17.9	
$32.0 \times 10^{-5}$ M	25.0	
$24.0 \times 10^{-5}$ M	35.7	
$16.0 \times 10^{-5}$ M	50.2	
$8.0 \times 10^{-5}$ M	70.8	

2. What would be the effect of the following on the Beer's Law plot?
- At the higher concentrations used to construct the plot, some solute molecules react with one another to form a species that doesn't absorb visible light.
  - A fingerprint is present on the solution cell where the light beam passes through.
  - Some solute molecules react with one another to form a species that doesn't absorb visible light.

## COLORIMETRY: A BEER'S LAW INVESTIGATION

### Tips and Traps

1. Page 3 of the experiment has a large blank space where students are to make notes on the development of Beer's Law as provided by the instructor. One suggestion for development follows.

A plot of  $\text{Log}(T)$  (y-axis) vs. concentration is linear, but the y intercept of the plot is not zero. A zero y intercept might be desired so that the measured function will be zero at zero concentration. For  $\text{Log}(\text{number})$  to be zero, number must equal exactly one. To obtain a number exactly equal to one the easiest thing to do is divide the measured number by itself. At zero concentration the measured current is the blank current  $I_0$ . Hence the solution to the intercept problem is to divide each current measurement by  $I_0$  and then take the log of this ratio, i.e.,  $\text{Log}(I/I_0)$ . **Have the students do this. NOW!** Use the **Add Formula** to calculate  $\text{Log}(\text{current}(I)/\text{measured value of } I_0)$ . For example, if  $I_0$  was measured to be 370.23 microamps, then the function for Column C would be  $\text{Log}(I/370.23)$ . "Click-drag" this formula to Column C.

Also note that the slope of the line is negative, not positive. A positive slope might be desired so that the measured function will increase as concentration increases. To change the slope from negative to positive is easy. Simply multiply the function by -1. **Have the students do this NOW** by adding a new formula. "Click-drag" this formula to Column D of the spreadsheet to be  $-\text{Log}(I/I_0)$  and to the **Y2** graph axis.

The plot of Column D vs. concentration is now linear with a zero intercept and a positive slope. The function  $-\text{Log}(I/I_0)$  is called the absorbance and Beer's Law has been established.

2. Since the concentration of the crystal violet solutions are so small, the concentration column should be formatted in Scientific Notation. This will allow proper printing of the data table from the Spreadsheet. To cause the spreadsheet to record the input in **scientific notation**, **right click** on that spreadsheet column and select **Column Properties**, and click on **Scientific Notation**.
3. Crystal violet stains, so appropriate care should be taken in its handling.
4. Be sure to caution the students to **PLACE THE CAP** before clicking **Add**.
5. Be sure the students mark the vial so it can be placed in the colorimeter chamber the same way each time.

## COLORIMETRY: A BEER'S LAW INVESTIGATION

### Suggested Answers for the Report Sheet

Dilution Table for Standards

	<u>#1</u>	<u>#2</u>	<u>#3</u>	<u>#4</u>	<u>#5</u>
<u>Desired CV</u>	<u>2.00 x 10<sup>-6</sup></u>	<u>4.00 x 10<sup>-6</sup></u>	<u>6.00 x 10<sup>-6</sup></u>	<u>8.00 x 10<sup>-6</sup></u>	<u>10.0 x 10<sup>-6</sup></u>
<u>ml CV</u>	<u>3.00</u>	<u>6.00</u>	<u>9.00</u>	<u>12.00</u>	<u>15.00</u>
<u>ml H<sub>2</sub>O</u>	<u>12.00</u>	<u>9.00</u>	<u>6.00</u>	<u>3.00</u>	<u>0.00</u>
<u>Total ml</u>	<u>15.00</u>	<u>15.00</u>	<u>15.00</u>	<u>15.00</u>	<u>15.00</u>

1. Record the number of your unknown \_\_\_\_\_

Record the absorbance of your unknown \_\_\_\_\_

Record the concentration of your unknown \_\_\_\_\_

2. Show your calculations for the concentration of your unknown.

$$A = \epsilon bc, \text{ Therefore, } c = A/(\epsilon b)$$

**Students will use their own values of molar absorptivity " $\epsilon$ " and absorbance "A." The value of b is 2.54 cm.**

3. As discussed by your instructor, Beer's Law is expressed by the equation  $A = \epsilon bc$ . In this equation, A is the absorbance of the solution, b is the cell path length (2.54 cm in this case), c is the molar concentration of the light-absorbing solute (CV), and  $\epsilon$  is a proportionality constant called the molar absorptivity.

a. In the space below, calculate the molar absorptivity for each of the five standard CV solutions and determine the average of the five values. Watch significant figures and be sure to include proper units. Note how well your values agree with one another. After all,  $\epsilon$  is supposed to be a constant!

**$\epsilon = A/bc$  where  $b = 2.54$  cm. Since "b" is a constant, and A is proportional to c, then  $\epsilon$  should be a constant for each solution. The values will vary as a function of the experimental conditions.**

b. Compare the numerical value of your average molar absorptivity with the slope of the best straight line given on your Beer's Law plot. That slope represents the "best" estimate of  $\epsilon$  from your experimental data.

\_\_\_\_\_ **The average of the five molar absorptivity values should be very close to the slope of the Beer's Law plot.**

## COLORIMETRY: A BEER'S LAW INVESTIGATION

### Suggested Answers to Questions/Problems

1. When iron (III) chloride is added to a solution of salicylic acid a purple solution forms. The photocell current read with the colorimeter for a series of standard solutions is given below. Convert these values to absorbencies. (The blank current is 100.0.)

FeCl <sub>3</sub>	Transmittance	Absorbance
40.0 x 10 <sup>-5</sup> M	17.9	0.747
32.0 x 10 <sup>-5</sup> M	25.0	0.602
24.0 x 10 <sup>-5</sup> M	35.7	0.447
16.0 x 10 <sup>-5</sup> M	50.2	0.299
8.0 x 10 <sup>-5</sup> M	70.8	0.150

2. What would be the effect of the following on the Beer's Law plot?
- a. At the higher concentrations used to construct the plot, some solute molecules react with one another to form a species that doesn't absorb visible light.

*This factor produces non-linear plot behavior.*

- b. A fingerprint is present on the solution cell where the light beam passes through.

*The Beer's Law plot line is displaced toward higher absorbance values.*

- c. Some solute molecules react with one another to form a species that doesn't absorb visible light.

*This factor produces non-linear plot behavior.*

# COLORIMETRY: A BEER'S LAW INVESTIGATION

## Sample Graphs




Fig. 1. Current versus Concentration.




Fig. 2. Current squared versus Concentration.




Fig. 3. 1/Current versus Concentration.




Fig. 4. Log(Current) versus Concentration




Fig. 5. -Log(Current) versus Concentration.




Fig. 6. -Log(Current/370.23) versus Concentration



## COLORIMETRY: A BEER'S LAW INVESTIGATION

### Laboratory Preparation (per student station)

#### Equipment

- **MicroLAB** colorimeter
- 6 colorimeter vials, 2.54 cm path length
- ring stand
- buret clamp
- buret with at least 15 ml capacity

#### Supplies

- KimWipes

#### Chemicals

Actual quantities needed are given below. A 50% excess is recommended.

- $1.00 \times 10^{-5}$  M crystal violet solution (45.00 ml)

#### Safety and Disposal

- crystal violet may cause skin and eye irritation and the NaOH solutions are corrosive. Make sure students wear goggles at all times.
- all resulting solutions may be flushed down the drain with plenty of water.
- have students wash hands with soap and water before leaving the lab.

