

What are computer requirements for MicroLab?



Minimum requirements:

The latest version of MicroLab software is compatible with 64 bit and 32 bit Windows 7 and Windows Vista as well as with Windows XP, Windows 2000, and Windows 98.

The hardware requirements are Pentium III or later (or analogous AMD processor) or Intel powered MAC with XP or Vista, 512 MB RAM, 16 MB of hard drive space, and a USB port to connect our FS522 lab interface. Our software runs just fine on netbooks with flash drives and Atom processors.

In all likelihood, the minimum hardware and operating system requirements of your institution or of your other software (such as Office™ or web browsers) will exceed those of MicroLab. Data acquisition and graphing is not a computer-intensive operation.



MicroLab software is compatible with 64 and 32-bit Windows 7, Windows Vista, Windows XP, and with Mac's running Intel emulators.
